



Fresh

Osteochondral Plugs

Clinical Overview

Fresh Osteochondral Plugs are pre-cut and viable osteochondral allografts for restoring small osteochondral defects. The hyaline cartilage and intact bony scaffold provide a durable and natural restoration of the damaged area.

Applications

- Small osteochondral defect repair

Why Use

- **Osteoconductive:** Created to precisely fit the recipient site when using the Arthrex® OATS® system
- **Compatible:** When used in conjunction with the Arthrex OATS system, the fresh plug will provide a line-to-line fit in the donor site. Graft can be easily downsized in diameter utilising the Arthrex OATS kits.
- **Patient-Friendly:** Eliminates donor site morbidity from allograft recovery. May decrease OR time due to decreased graft prep time. Less OR time can mean less time under anesthesia and less tourniquet time.^{1, 2}
- **Viable chondrocytes:** Validated aseptic processing and never frozen to ensure chondrocyte viability

References

1. Boniello M, Robinson SP, Bonner KF. Osteochondral Allograft Transplantation. Chapter 31: Insall and Scott Surgery of the Knee Textbook. 6th ed. Edited by W. Norman Scott; Elsevier, Churchill, & Livingstone. Philadelphia, PA. 2017: 442-453.
2. Poehling G, Harner C, Verdonk R. Symposium: Allografts in sports medicine. Presented at the 7th European Federation of National Associations of Orthopaedics and Traumatology Congress. June 4-7, 2005. Lisbon.

Speak to your local Business Development Manager,
for further information using the details below:


T: 01443 719 555

E: customerservice@hospitalinnovations.co.uk

W: www.hospitalinnovations.com

Hospital Innovations Limited
Concept House
Talbot Green Business Park
Pontyclun
CF72 9FG

HOSPITAL
INNOVATIONS



Fresh Osteochondral Plugs with viable cartilage

Code	Diameter (D)
RFP10	10 mm
RFP16	16 mm

Set length of 13± 1 mm • Refrigeration storage (1-10 °C) • 45 day shelf life

Speak to your local Business Development Manager,
for further information using the details below:

T: 01443 719 555

E: customerservice@hospitalinnovations.co.uk

W: www.hospitalinnovations.com

Hospital Innovations Limited
Concept House
Talbot Green Business Park
Pontyclun
CF72 9FG

HOSPITAL
INNOVATIONS