


# Costal Cartilage

## A natural, biological solution



Costal Cartilage is sterile (SAL  $10^{-6}$ ) hyaline cartilage from the rib. Cleaned and sterilised using the patented Allowash XG® technology, this graft's biomechanical and biochemical properties are maintained, making it ideal for nasal and ear reconstruction.

Product Name	<ul style="list-style-type: none"><li>• Costal Cartilage</li></ul>
Clinical Applications	<ul style="list-style-type: none"><li>• Nasal Reconstruction</li><li>• Ear Reconstruction</li></ul>
Features And Benefits	<ul style="list-style-type: none"><li>• <b>Sterile:</b> Sterilised using proprietary and patented Allowash XG technology. This technology provides a sterility assurance level of <math>10^{-6}</math>, without compromising the graft's inherent properties</li><li>• <b>Biocompatible:</b> Costal Cartilage is composed of 100% hyaline cartilage that is remodeled naturally during the healing process</li><li>• <b>Supportive:</b> Provides biomechanical strength requirements for its intended use</li><li>• <b>Convenient:</b> Allows the surgeon to select the most appropriate graft for clinical procedure</li><li>• <b>Malleable:</b> Easy to trim and shape to perfectly fit surgical need</li></ul>



## Costal Cartilage

Size	Preservation Method	Order No.
L: 30 - 59 mm; D: $\geq 7$ mm; W: $\geq 7$ mm	Frozen	CCART30
L: 60 - 90 mm; D: $\geq 7$ mm; W: $\geq 7$ mm	Frozen	CCART60

Speak to your local Sales Specialist for further information or contact us using the details below:

**T:** 01443 719 555

**E:** [info@hospitalinnovations.co.uk](mailto:info@hospitalinnovations.co.uk)

[www.hospitalinnovations.com](http://www.hospitalinnovations.com)

Hospital Innovations Limited  
 Concept House  
 Talbot Green Business Park  
 Pontyclun  
 CF72 9FG

**HOSPITAL**  
**INNOVAT10NS**  
 Celebrating a decade of innovation