

# Allograft CANNULATED CLOWARD DOWELS

Maxxeus Allograft Cannulated Cloward Dowels offers surgeons a rapid and effective solution for filling bone tunnels during ACL/PCL revisions. The use of Cannulated Cloward Dowels provides immediate structural architecture for stability. Surgeon has the choice to rehydrate the dowels in the solution of their choice for the best incorporation of the platform.



- Maxxeus Allograft Cannulated Cloward Dowels have been processed to remove blood and bone marrow components.
- The dowels are pre-shaped with a tapered tip and cannulated for easier graft insertion.
- Maxxeus Allograft Cannulated Cloward Dowels can be stored in their packaged state for up to 5 years.
- Maxxeus Allograft Cannulated Cloward Dowels are terminally sterilized to achieve a Sterility Assurance Level (SAL) of 10<sup>-6</sup> within the final packaging.

Allograft	Size	Graft Code	Preservation
Cannulated Cloward Dowel, 8 mm	8 x 26mm+	3617-12	Freeze Dried
Cannulated Cloward Dowel, 9 mm	9 x 26mm+	3618-12	Freeze Dried
Cannulated Cloward Dowel, 10 mm	10 x 26mm+	3626-12	Freeze Dried
Cannulated Cloward Dowel, 11 mm	11 x 26mm+	3619-12	Freeze Dried
Cannulated Cloward Dowel, 12 mm	12 x 26mm+	3620-12	Freeze Dried
Cannulated Cloward Dowel, 14 mm	14 x 26mm+	3621-12	Freeze Dried
Cannulated Cloward Dowel, 16 mm	16 x 26mm+	3622-12	Freeze Dried
Cannulated Cloward Dowel, 18 mm	18 x 26mm+	3623-12	Freeze Dried

Speak to your local Business Development Manager, for further information using the details below:

T: 01443 719 555

E: [customerservice@hospitalinnovations.co.uk](mailto:customerservice@hospitalinnovations.co.uk)

W: [www.hospitalinnovations.com](http://www.hospitalinnovations.com)

Hospital Innovations Limited  
Concept House  
Talbot Green Business Park  
Pontyclun  
CF72 9FG