

ReadiGRAFT® Cancellous Chips



ReadiGRAFT® Cancellous Chips are a natural osteoconductive scaffold. Chips are available in 4 grind sizes and a variety of volumes to provide maximum surgical flexibility.

| | |
|------------------------------|--|
| <p>Clinical Applications</p> | <p>Any surgical application that requires bone void filler</p> |
| <p>Features And Benefits</p> | <p>Natural Osteoconductive Scaffold: The natural interconnected porous architecture provides an optimal scaffold for cellular attachment, proliferation, and vascular in-growth.</p> <p>Hydrophilic: Absorbs and retains bioactive fluids like blood, Platelet Rich Plasma (PRP), and Bone Marrow Aspirate (BMA).</p> <p>Sterile: The implant is sterilised using proprietary and patented Allowash XG® technology. This technology provides a sterility assurance level of 10⁻⁶, without compromising the implant's inherent osteoconductivity or osteoinductive potential.</p> |
| <p>Alternative Products</p> | <p>Cancellous Cubes Cortical Cancellous Chips Femoral Head (Grinder)</p> |

ReadiGRAFT®

Cancellous Chips

Freeze-dried Chips

| Item No. | Item Description |
|-------------|---|
| CAN15 14BP | Freeze-dried Crushed Cancellous Chips: 1-4mm, 15cc |
| CAN30 14BP | Freeze-dried Crushed Cancellous Chips: 1-4mm, 30cc |
| CAN60 14BP | Freeze-dried Crushed Cancellous Chips: 1-4mm, 60cc |
| CAN5 | Freeze-dried Crushed Cancellous Chips: 1-8mm, 5cc |
| CAN15 | Freeze-dried Crushed Cancellous Chips: 1-8mm, 15cc |
| CAN30 | Freeze-dried Crushed Cancellous Chips: 1-8mm, 30cc |
| CAN60 | Freeze-dried Crushed Cancellous Chips: 1-8mm, 60cc |
| CAN15 410BP | Freeze-dried Crushed Cancellous Chips: 4-10mm, 15cc |
| CAN30 410BP | Freeze-dried Crushed Cancellous Chips: 4-10mm, 30cc |
| CAN60 410BP | Freeze-dried Crushed Cancellous Chips: 4-10mm, 60cc |
| CAN90 410BP | Freeze-dried Crushed Cancellous Chips: 4-10mm, 90cc |

Frozen Chips

| Item No. | Item Description |
|-----------|--|
| FCAN30 14 | Frozen Crushed Cancellous Chips: 1-4mm, 30cc |
| FCAN60 14 | Frozen Crushed Cancellous Chips: 1-4mm, 60cc |
| FCAN30 | Frozen Crushed Cancellous Chips: 1-8mm, 30cc |

Speak to your local Business Development Manager,
for further information using the details below:

T: 01443 719 555
E: customerservice@hospitalinnovations.co.uk
W: www.hospitalinnovations.com

Hospital Innovations Limited
Concept House
Talbot Green Business Park
Pontyclun
CF72 9FG