



# OraGRAFT®

## Ribs

**Clinical Overview** The OraGraft allograft line is comprised of hard and soft allograft implants generally used for sinus augmentation, craniofacial reconstruction and correction of periodontal and ridge defects

- Applications**
- Ridge Maintenance
  - Craniofacial Reconstruction

- Features & Benefits**
- **Highest Safety Level:** Sterile with an SAL of  $10^{-6}$ .
  - **Convenient Storage:** Ambient temperature with a 5 year shelf life.
  - **Ease of Use:** Can be affixed to surgical site using standard bone screws.
  - **Versatile:** Solid cortical shell with cancellous contained within.
  - **Match to Defect:** Easy to trim, can be shaped using standard instrumentation.



## OraGraft Ribs

Ambient Temperature\*/5 year shelf life

Sizes	Order Code
60 - 115 mm	RIB

\*While ambient room temperature has not been defined by regulatory bodies, LifeNet Health would recommend storage at 2°C to 37°C with excursions of less than 24 hours up to 40°C. If an excursion outside this range occurs, please contact LifeNet Health.

Instructions for use available at [LifeNetHealth.org/IFU](http://LifeNetHealth.org/IFU)

Speak to your local Business Development Manager for further information or contact us using the details below:

**T:** 01443 719 555

**E:** [info@hospitalinnovations.co.uk](mailto:info@hospitalinnovations.co.uk)

[www.hospitalinnovations.com](http://www.hospitalinnovations.com)

Hospital Innovations Limited  
Concept House  
Talbot Green Business Park  
Pontyclun  
CF72 9FG

HOSPITAL  
INNOVATIONS