


NanoFx[®]

A-Curve

Smaller. Deeper. Better.

“The overriding finding of the present study is that small subchondral drill holes reflecting the physiological subchondral trabecular distance significantly improve osteochondral repair.”*

Eldracher M, et. al. Small subchondral drill holes improve marrow stimulation of articular cartilage defects. AJSM 2014 Nov;42(11):2741-50.



- Specially designed curve for hard to access ankle lesions
- Reduces surface skyving for off-axis channel preparation
- Improves bone marrow access for cartilage repair procedures



Thumble



Handle Instrument




PleuriStik™ Guide Wire

Item No.	Item Description
FURS-010R	NanoFx [®] Thumb Tab Accessory
5500-4010	NanoFx [®] A Curve Handle
5500-1020	NanoFx [®] Handle 15
FURS-2101	NanoFx [®] PleuriStik 15 - Box of 5
FURS-4101	NanoFx [®] PleuriStik A - Box of 5
FURS-2101/1	NanoFx [®] PleuriStik 15 - Single
FURS-4101/1	NanoFx [®] PleuriStik - Single


Smaller

“Significant enhancements were observed at the level of individual parameters and of overall histological articular cartilage repair, together with improved immunoreactivity to type II and type I collagen of the cartilaginous repair tissue. Second, the microarchitecture of both the subchondral bone plate and the subarticular spongiosa was better reconstituted.”*

Eldracher M, Orth P, Cucchiari M, Pape D, Madry H. Small subchondral drill holes improve marrow stimulation of articular cartilage defects. Am J Sports Med. 2014 Nov;42(11):2741-50.




DEEPER



STEP 1

Assemble & orient on lesion



STEP 2


Gently tap to depth stop



STEP 3

Push thumb to extract pleuristik. Repeat.

BETTER



Microfracture

1mm K-Wire

NanoFracture

Figures: open trabecular channels; closed trabecular channels, microCT comparison: Axial (top), Sagittal (bottom).

Speak to your local Business Development Manager, for further information using the details below:

T: 01443 719 555

E: customerservice@hospitalinnovations.co.uk

W: www.hospitalinnovations.com

Hospital Innovations Limited
 Concept House
 Talbot Green Business Park
 Pontyclun
 CF72 9FG